

How to Run Patches in WinRecs 3



MED2020 Health Care Software Inc.

MED2020 Health Software Inc. (MED2020) is the leading provider of modular health information management solutions for the Canadian health care industry and is headquartered in Ottawa, Ontario. Solutions are offered to assist in capturing, reporting and analyzing health data to support enhanced information sharing, encourage informed decision making and streamline facility operations. MED2020's flagship product, WinRecs, is the foundation for a complementary suite of modules that have also been designed to assist health information management departments with their operational needs.

Contact Details

MED2020 is open Monday to Friday between 07:00 and 19:00 (Eastern Time), except Canadian federal holidays and Ontario provincial holidays.

Head Office

MED2020 Health Care Software Incorporated
4471 Innes Road, Suite #200
Ottawa, Ontario
K4A 3J7

Tel: (613) 830-3761


Fax: (613) 830-2410

Toll Free: (800) 461-2020


E-mail: support@med2020.ca

Web: <http://www.med2020.ca>

Step 1 : Launch SQL Management Studio





Step 2: Connect to the WinRecs database server



Step 3: Select the WinRecs patch to run


The file will have an extension of .SQL






Step 4: Select the database to update


Choose the applicable WinRecs database from the drop down list



Step 5: Click execute



The script will run and the results will display in the message window
If the results state that the 'query executed successfully' then the process is complete



If the results state 'Query completed with errors' contact 1-800-461-2020 or support@med2020.ca for further instruction

