





WinRecs 3 Version 2024.03 Interface Services Update

In comparison to previous WinRecs version updates, this update requires additional steps outlined below. The changes are used for the configuration of the following MED2020 services.




The steps require downloading files and placing them in the following folders on the server where your MED2020 interface services are installed.



MED2020 Listener Service

For clients that have a **HL7 interface** and the 'MED2020 HL7 Listener' and the 'MED2020 HL7 Queue' service installed, first download this file, unzip it and paste the folder into the 'HL7 Listener' like this screenshot.


https://files.med2020.ca/Downloads/Winrecs3/Winrecs_Releases/2024.03/Interfaces/HL7%20Listener/HL7%20Listener%203.24.3.12.zip



MED2020 Queue Service

And then download this file, unzip it and paste the folder into the 'HL7 Queue' Folder

https://files.med2020.ca/Downloads/Winrecs3/Winrecs_Releases/2024.03/Interfaces/HL7%20Queue/HL7%20Queue%203.24.3.12.zip




The screenshot shows a Windows File Explorer window with the address bar set to 'This PC > Local Disk (C:) > ProgramData > MED2020 > HL7 Queue >'. The main pane displays a table of folders:

Name	Date modified	Type	Size
3.24.2.28	6/11/2024 3:00 PM	File folder	
3.24.2.95	7/19/2024 12:32 PM	File folder	
3.24.3.12	10/10/2024 10:08 AM	File folder	

MED2020 Autobatch Service

For clients that have a **scheduled batch interface** and the 'MED2020 Auto Batch Interface' service installed, download this file, unzip it and paste the folder into the 'Auto batch folder' like this screenshot.

https://files.med2020.ca/Downloads/Winrecs3/Winrecs_Releases/2024.03/Interfaces/Autobatch/Autobatch%203.24.2.4.zip



The screenshot shows a Windows File Explorer window with the address bar set to 'This PC > Local Disk (C:) > ProgramData > MED2020 > Auto Batch >'. The main pane displays a table of folders:

Name	Date modified
3.24.3.12	10/10/2024 12:43 PM

For assistance contact support@med2020.ca